

Instruction for use

myQuickmat Classico kit

 polydentia
swiss manufacture 

 sectional
matrix
systems

Medical Devices Kit for professional use only - REV 02 06.2021

instruction for use

A - Description

myQuickmat Classico Kit: Complete sectional matrix system developed to offer gentle tooth separation in order to create contact points in the restoration of class II cavities. This kit features Polydentia's myRing Classico, a classic sectional matrix ring made of nickel titanium. It is reinforced to provide an optimal and gentle separation force in combination with the included Interdental Maple Wood Wedges. For the restoration of very extensive cavities, the kit also includes the newly re-designed Delta Tubes. These anatomically shaped silicone tubes ensure optimal matrix retention and adaptation to the tooth's anatomy. LumiContrast sectional matrices complete the kit's assortment. These dark-blue matrices significantly increase contrast and reduce glare effect, especially while using loupes or microscopes.

B - Composition

myRing Classico: NiTi Ring with thermoplastic overmolding
Delta Tubes: silicone
LumiContrast matrices: Stainless steel
Quickmat matrices: Stainless steel
Wood wedge: Maple wood

C - Indications

myQuickmat Classico Kit is indicated for gentle tooth separation in order to create correct contact points in the restoration of class II cavities.

D - Contraindications

None Known.

E - Side effects

There are no known side effects. Promptly report any serious incident by contacting your distributor or Polydentia.

F - Interactions with other medical devices or products

There are no known interactions.

G - Information for the patient

No information for the patient related to the myQuickmat Classico Kit and its component is necessary.

H - Shelf life

Please refer to the expiration date on the packaging

I - Preservation and storage

Store in a clean, dry environment away from moisture or excessive humidity.

L - Warnings

Due to the presence of small and sharp components, the use of a rubber dam is mandatory to avoid accidental ingestion or aspiration. Always ensure the rings using dental floss. In case of accidental inspiration of a myQuickmat Classico Kit component (matrix, retaining ring and its extremities, wedge) seek immediate medical attention. In case of ingestion of a myQuickmat Classico Kit component, it is strongly recommended seeking for medical assistance and locate the ingested components. If illness develops after ingestion, consult a doctor immediately. All Polydentia reusable Products are supplied non-sterile and must be appropriately disinfected, cleaned and sterilized before the first use. Please refer to the Cleaning, disinfection and Sterilization Instructions at www.polydentia.com.

M - Application

See the following pictograms.

N - Pack sizes available

See "Reorder now" section.

O - Disposal

Because of their composition, Matrices and wedges can not be sterilized and are therefore single use devices. Do not reuse them. Properly dispose devices contaminated with potentially infectious substances of human origin, in accordance with local regulations.




P - Ring reconditioning

For a proper, effective and safe reconditioning of the myRing Classico and their extremities read the Cleaning, Disinfection and Sterilization Instructions at www.polydentia.com.

QR instructions
Cleaning and
Sterilization
Guide or go to
www.polydentia.ch

application

1. Prepare the cavity

Remove the carious lesion and prepare tooth with a conservative Class II preparation. Even if the contact was not totally broken the matrix can be placed.

2. Place the matrix

Select an appropriate matrix band according to the tooth morphology and restoration type. Place the contoured sectional matrix with the concave edge towards the occlusal margin using your fingers or with tweezers.

3. Wedging

Select the appropriate sized wooden wedge to close the gingival margin, and insert it firmly. When needed maintain a finger on the matrix band to avoid dislodging during wedging.

4. Place myRing Classico

Expand myRing Classico using Polydentia's myForceps or standard rubber dam forceps. Place the ring carefully and securely. Please refer to the Tips & Tricks section to see the different ways to position myRing Classico and ensure best teeth separation as well as matrix adaptation.

5. Burnish the matrix

For an excellent contact burnish the band against the adjacent tooth using a ball burnisher or tweezer.

6. Restore the cavity

Restore the cavity according to the filling material manufacturer's instructions.

7. Remove the sectional matrix system

Remove myRing Classico, the wedge and the matrix. Our suggestion for an easy matrix removal is to use the special Polydentia myQuickmatrix Forceps.

8. Finishing and polishing

With thin burs remove any small excess of material then polish the restoration according to the finishing and polishing instruments manufacturer's recommendations.


tips & tricks using myQuickmat Classico kit

Wedging


According to the clinical situation, the wooden wedge can be broken in order to better accommodate myRing Classico with Delta Tubes. This trick is very helpful in case of small teeth.

Delta Tubes placement options

Delta Tubes can be oriented according to the clinical situation to obtain the best matrix retention and adaptation to the tooth's anatomy. The most common cases are explained as follows (Case A, B). Different orientation possibilities can be combined to best adapt myRing Classico to each specific clinical situation. The use of Delta Tubes is not mandatory, myRing Classico can be placed also without Delta Tubes (Case C, D).

Case A, small cavities:

Delta Tubes edges allow a perfect adaptation between the teeth

Case B, broad cavities:

Delta Tubes flat surfaces provide support on the adjacent tooth, no matrix deformation.

Case C and D, myRing Classico without Delta Tubes


myRing Classico tines placed near to the wood wedge

This placement technique is suitable for cases of a wide space in the interproximal area.

myRing Classico tines placed over the wood wedge

This placement technique is suitable if the interproximal space is tight or when the cavity is deep.

MOD and multiple tooth restorations:

MyQuickmat Classico rings are designed to be stacked over one another to perform a MOD and multiple tooth restorations.

1. Prepare the cavities as described in the instructions for use, step 1.
2. Place matrices and wedges on both mesial and distal surfaces of the tooth being restored
3. Position myRing Classico
4. Complete the restoration as described in the instructions for use, steps 5 to 8.

symbols

REF

Item number

LOT

Batch code

Use by

Consult instructions for use

Do not reuse

Caution

Rx Only

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner

Manufacturer

CE mark

Ultrasonic bath

Washer-disinfector for thermal disinfection

Sterilizable in a steam sterilizer (autoclave) at the temperature specified

MD

Medical Device

EC REP

European Authorised Representative

Polydentia SA
via Cantonale 47
6805 Mezzovico-Vira
Switzerland

phone: +41 91 946 29 48
fax: +41 91 946 32 03

info@polydentia.ch
www.polydentia.ch

reprocessing of Polydentia reusable products

Disinfection,
Cleaning &
Sterilisation
instructions

Index

<u>EN - English</u>	2
<u>DE - Deutsch</u>	6
<u>FR - Français</u>	11
<u>IT - Italiano</u>	16
<u>ES - Español</u>	21

Reprocessing of Polydentia reusable Products

Disinfection, Cleaning & Sterilisation instructions

Basic Principles:

- Polydentia Products are supplied non-sterile and must be appropriately disinfected, cleaned and sterilized before the first use.
- Polydentia Products must be disinfected, cleaned and sterilized prior to each use.
- Effective cleaning is an indispensable requirement for proper Polydentia Products sterilization.
- Consider the legal provisions valid for your country as well as to the hygienic instructions of the Dental Practice.
- The user is responsible for the sterility of Polydentia Products. Therefore, please ensure that only validated procedures are used for disinfection, cleaning and sterilization. The sterilization equipment must also be maintained and checked regularly, as well as the validated parameters applied to each disinfection, cleaning and sterilization cycle.
- Materials may alter over time. Sterilization or exposure to chemicals can accelerate this deterioration. Always check your Polydentia Product for wear and damage before use.
- Avoid contact of dissimilar metals at any point during the disinfection, cleaning and sterilization process.

Protection of Staff Members:

Used and contaminated Polydentia Products must be handled with appropriate Personal Protective Equipment.

1. Pre-treatment

Before processing Polydentia Products remove coarse impurities on the products immediately after application.

Do not use a metal brush. A metal brush will cause scratches or damages on the products surface.

- In case of using a thermal disinfectant the steps in this guidance (see below) “2. Disinfection” and “3. Rinsing” are unnecessary, provided that the dirty Polydentia Products have not been stored dry longer than 2 hours.

2. Disifencion

- Do not leave dirty Polydentia Products dry – put them into a disinfecting liquid as soon as possible.
- We suggest to use an enzymatic pre-cleaner disinfectant such as, but not limited to:
 - Dürr Dental ID215
 - Or a pre-cleaning disinfectants based on Quaternary Ammonium compounds.
- Follow carefully the instructions for use given by the manufacturer of the selected disinfectant. Please note that too strong concentration of the disinfectant or too long exposure time may damage the products.
- Change the disinfectant regularly according to the instructions.

3. Rinsing

- Immediately after the disinfection, Polydentia Products must be rinsed carefully under running lukewarm water.

4. Cleaning

Always clean Polydentia Products, disinfection and rinsing are not sufficient.

Cleaning alternatives:

Manual (A), Ultrasonic bath (B) Thermal Disinfector (C)

If possible, an automatic procedure in a thermal disinfector should be used for cleaning of Polydentia Products. A manual procedure, even in case of application of an ultrasonic bath, should only be used if an automatic procedure is not available or if such a method is not compatible with specific materials; in this case, the significantly lower efficiency of a manual procedure must be considered.

A. Manual Cleaning

- Completely disassemble Polydentia Products if applicable
- We suggest to use a cleaning solution such as, but not limited to:
 - Dürr Dental ID212Or cleaning solutions based on Quaternary Ammonium compounds. Check the pH-value of the cleaning solutions; it must be as neutral as possible.
- Follow carefully the instructions for use given by the manufacturer of the selected cleaning solutions. Please note that too strong concentration or too long exposure time may damage Polydentia Products.
- Change the cleaning solutions regularly according to the instructions.
- Always rinse carefully Polydentia Products under running water after cleaning.

B. Ultrasonic bath

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Completely disassemble Polydentia Products if applicable.
- Follow carefully the instructions for use of the ultrasonic bath device.
- Service and clean the device frequently, according to the instructions for use.
- Ensure that maximum load is not exceeded.
- We recommended to use a cleaning solutions such as, but not limited to:
 - Dürr Dental ID212
 - Or cleaning solutions based on Quaternary Ammonium compounds. Check the pH-value of the cleaning solutions; it must be as neutral as possible.
- Follow carefully the instructions for use given by the manufacturer of the selected cleaning solutions. Please note that too strong concentration or too long exposure time may damage Polydentia Products.
- Change the cleaning solutions regularly according to the instructions.
- Remove the devices from the ultrasonic bath immediately after cleaning.
- Always rinse carefully Polydentia Products under running water after cleaning.

C. Thermal Disinfector

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Completely disassemble Polydentia Products if applicable.
- Follow carefully the instructions for use given by the manufacturer of the thermal disinfector and of the selected program; and prefer cleaning agents with corrosion protector.
- Inspect the device regularly according to the instructions for use. Service and clean the device frequently, according to the instructions for use.
- Ensure that maximum load is not exceeded.
- Ensure that the selected washing programme and cleaning agents are effective enough.
- Remove Polydentia Products from the thermal disinfector immediately after cleaning.

5. Drying

Before sterilizing in the autoclave, Polydentia Products must be thoroughly dried. The normal tap water may contain different chemical ingredients which may cause chemical reactions when autoclaving.

6. Checking

After cleaning, check that Polydentia Products are clean by rubbing the surface with an eraser or finger.

WARNING! It is extremely important to check that the products are truly clean and dry before the sterilizing procedure.

7. Sterilization

Sterilize Polydentia Products **only** in a steam autoclave using distilled water and at the temperature specified on the relevant symbol. **Do not use chemical or cold sterilization.**

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Follow carefully the instructions for use of the autoclave.
- Inspect the device regularly according to the instructions for use. Service and clean the device frequently, according to the instructions for use .
- Ensure that maximum load is not exceeded.
- Place all the medical device components in an autoclave tray, or in a sterilization disposable pouch (recommended).
- Plastic parts must not touch the autoclave walls, since the temperature there might be higher.
- We recommend to use always the following cycle (Prion Cycle):
Temperature of sterilisation: 134°C
Pressure of sterilisation: 2 bar
Time of sterilisation: 20.5 min
Drying time: 20 min
- Remove Polydentia Products from the autoclave immediately after sterilisation.

8. Storage

- Store sterilized product in a dry and dust free location.
- Keep product sealed in the sterilization pouch until ready to use.

DISCLAIMER:

This high quality product is manufactured by Polydentia SA, in Mezzovico, Switzerland. The content of this packaging has been inspected, tested and approved by our quality officers.

Materials used during our manufacturing process are designed to provide durability and performance. Polydentia constantly looks for the highest quality plastics for its product; nevertheless please take into account that plastics may alter over time. Sterilization or exposure to chemicals can accelerate this deterioration.

For the longest lifetime of your premium quality product always clean and sterilize it following these instructions.

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