Instruction for use

unica introkit anterior

developed with
STYLE ITALIANO™

the one solution
for anterior
restorations

 polydentia
swiss manufacture 

 matrix
systems

instruction for use

A - Description

Unica anterior and Unica minideep anterior, thanks to their contoured shape, adapt correctly to the different morphologies of anterior teeth and are intended to restore interproximal and cervical margins at once, even in presence of rubber-dam or gingival retraction cords reducing chair-time significantly. Furthermore, Unica anterior and Unica minideep anterior matrices, once positioned, allow to easily visualize the final shape of the restoration.

Unica anterior & Unica minideep anterior allow to:

- Restore both proximal and cervical margins at once
- Anatomically restore the proximal margins thanks to its contoured profile
- Easily manage the cervical area

B - Composition

Unica anterior & Unica minideep anterior: stainless steel

myWedge (XS, S, M): thermoplastic

myQuickmatrix Forceps: Stainless steel

C - Indications

The matricing solutions of Unica anterior introkit (Unica anterior matrix and Unica minideep anterior matrix) are indicated for:

- Class III cavity restorations
- Class IV cavity restorations
- Class V cavity restorations
- Direct stratification composite veneers
- shape modifications

Unica anterior is indicated for:

- maxillary central incisors

Unica minideep anterior is indicated for:

- maxillary lateral incisors
- mandibular central and lateral incisors
- conoid teeth
- peg laterals
- triangular teeth
- teeth with narrow cervical diameter.

A liquid dam can be used to secure Unica anterior and Unica minideep anterior in place in case of wide preparations.

myQuickmatrix Forceps ensures a firm grip of the matrix when removing it after the restoration.

D - Contraindications

No contraindications are known.

E - Side effects

No side effects known

F - Interactions with other medical devices or products

There are no known interactions.

G - Information for the patient

No information for the patient related to the unica anterior Introkit and its components is necessary.

H - Shelf life

Store in a clean, dry environment away from moisture or excessive humidity.

I - Preservation and storage

Due to the presence of small and sharp components, the use of a rubber dam is mandatory to avoid accidental ingestion or aspiration.

In case of accidental inspiration of a component of Unica anterior Introkit (matrix or wedge) seek immediate medical attention.

In case of ingestion of a component of Unica anterior Introkit component, it is strongly recommended seeking for medical assistance and locate the ingested components. If illness develops after ingestion, consult a doctor immediately.

myQuickmatrix forceps is supplied non-sterile and must be appropriately disinfected, cleaned and sterilized before the first use. Please refer to the Cleaning, disinfection and Sterilization Instructions at polydentia.ch

L - Application

See the following pictograms.

M - Pack sizes available

See "Order now" section.

N - Disposal

Because of their composition, matrices and wedges can not be sterilized and are therefore single use devices. Do not reuse them. Properly dispose devices contaminated with potentially infectious substances of human origin, in accordance with local regulations.

QR instructions
Cleaning and
Sterilization
Guide or go to
polydentia.ch

application

The following pictures show how to position and fix Unica anterior & Unica minideep anterior

select unica matrices

Choose Unica according to its shape and the tooth to be restored

Select *Unica anterior* for:

- maxillary central incisors
- maxillary lateral incisors
- canine teeth

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- maxillary and mandibular laterals
- mandibular central incisors
- conoid teeth
- triangular teeth
- peg laterals
- teeth with narrow cervical diameter

how to position Unica matrices

1 Use the placement wings to insert Unica proximally. In case of tight contact surfaces, insert one wing at a time.

2 Gently pull the matrix palatally by its placement wings

3 Push the cervical connection into the gingival sulcus

4 Insert myWedge to fix the matrix in position

tips and tricks

Secure Unica matrices in place with a liquid dam, in case a wedge risks displacing or deforming the matrix.

Yes

No

symbols

REF	Item number
LOT	Batch code
 i	Consult instructions for use
 ⊗	Do not reuse
 ■	Manufacturer
 CE	CE mark
 ⌚	Ultrasonic bath
 ☒	Washer-disinfector for thermal disinfection
 134°C	Sterilizable in a steam sterilizer (autoclave) at the temperature specified
 MD	Medical Device
EC REP	European Authorized Representative

 Polydentia SA
via Cantonale 47
6805 Mezzovico-Vira
Switzerland

phone: +41 91 946 29 48
fax: +41 91 946 32 03

info@polydentia.ch
polydentia.ch

 **polydentia**
swiss manufacture 

Unica Introkit anterior

REF 6900

20 pcs Unica anterior
 20 pcs Unica minideep anterior
 40 pcs myWedge XS
 40 pcs myWedge S
 40 pcs myWedge M
 1 pc myQuickmatrix Forceps

Use only with rubber dam - For professional use only

Refills

	Quantity	REF	CE	MD	⊗
	12	6912	CE	MD	⊗
	50	6950	CE	MD	⊗
	12	6911	CE	MD	⊗
	50	6951	CE	MD	⊗
	100	7101	CE	MD	⊗
	100	7102	CE	MD	⊗
	100	7103	CE	MD	⊗
	1	5703	CE	MD	⊗

Polydentia SA - via Cantonale 47
 6805 Mezzovico-Vira, Switzerland
 Made in Switzerland

phone: +41 91 946 29 48 - fax: +41 91 946 32 03
info@polydentia.ch - polydentia.ch
[f](#) [g](#) [in](#)

Information guide

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developed with
STYLE ITALIANO™

more information

REV 04.08.2021

All individual components are CE marked

matrix selection

Choose Unica according to its shape and the tooth to be restored

Select *Unica anterior* for:

- maxillary central incisors
- maxillary lateral incisors
- canine teeth

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- maxillary and mandibular laterals
- mandibular central incisors
- conoid teeth
- triangular teeth
- peg laterals
- teeth with narrow cervical diameter

clinical applications

Class III cavity restorations

Clinical image: courtesy of Dr Chiodera

Class IV cavity restorations

Clinical image: courtesy of Dr Jorquera

Class V cavity restorations

Clinical image: courtesy of Dr Chiodera

Shape modifications

Clinical image: courtesy of Dr Shaalan

Direct stratification composite veneers

Clinical image: courtesy of Dr Manauta

how to position Unica

1 Use the placement wings to insert Unica proximally. In case of tight contact surfaces, insert one wing at a time

2 Gently pull the matrix palatally by its placement wings

3 Push the cervical connection into the gingival sulcus

4 Insert a myWedge to fix the matrix in position

tips and tricks

Secure Unica matrices in place with a liquid dam in case a wedge risks displacing or deforming the matrix.

Yes

No

reprocessing of Polydentia reusable products

Disinfection,
Cleaning &
Sterilisation
instructions

Index

<u>EN - English</u>	2
<u>DE - Deutsch</u>	6
<u>FR - Français</u>	11
<u>IT - Italiano</u>	16
<u>ES - Español</u>	21

Reprocessing of Polydentia reusable Products

Disinfection, Cleaning & Sterilisation instructions

Basic Principles:

- Polydentia Products are supplied non-sterile and must be appropriately disinfected, cleaned and sterilized before the first use.
- Polydentia Products must be disinfected, cleaned and sterilized prior to each use.
- Effective cleaning is an indispensable requirement for proper Polydentia Products sterilization.
- Consider the legal provisions valid for your country as well as to the hygienic instructions of the Dental Practice.
- The user is responsible for the sterility of Polydentia Products. Therefore, please ensure that only validated procedures are used for disinfection, cleaning and sterilization. The sterilization equipment must also be maintained and checked regularly, as well as the validated parameters applied to each disinfection, cleaning and sterilization cycle.
- Materials may alter over time. Sterilization or exposure to chemicals can accelerate this deterioration. Always check your Polydentia Product for wear and damage before use.
- Avoid contact of dissimilar metals at any point during the disinfection, cleaning and sterilization process.

Protection of Staff Members:

Used and contaminated Polydentia Products must be handled with appropriate Personal Protective Equipment.

1. Pre-treatment

Before processing Polydentia Products remove coarse impurities on the products immediately after application.

Do not use a metal brush. A metal brush will cause scratches or damages on the products surface.

- In case of using a thermal disinfectant the steps in this guidance (see below) “2. Disinfection” and “3. Rinsing” are unnecessary, provided that the dirty Polydentia Products have not been stored dry longer than 2 hours.

2. Disifencion

- Do not leave dirty Polydentia Products dry – put them into a disinfecting liquid as soon as possible.
- We suggest to use an enzymatic pre-cleaner disinfectant such as, but not limited to:
 - Dürr Dental ID215
 - Or a pre-cleaning disinfectants based on Quaternary Ammonium compounds.
- Follow carefully the instructions for use given by the manufacturer of the selected disinfectant. Please note that too strong concentration of the disinfectant or too long exposure time may damage the products.
- Change the disinfectant regularly according to the instructions.

3. Rinsing

- Immediately after the disinfection, Polydentia Products must be rinsed carefully under running lukewarm water.

4. Cleaning

Always clean Polydentia Products, disinfection and rinsing are not sufficient.

Cleaning alternatives:

Manual (A), Ultrasonic bath (B) Thermal Disinfector (C)

If possible, an automatic procedure in a thermal disinfector should be used for cleaning of Polydentia Products. A manual procedure, even in case of application of an ultrasonic bath, should only be used if an automatic procedure is not available or if such a method is not compatible with specific materials; in this case, the significantly lower efficiency of a manual procedure must be considered.

A. Manual Cleaning

- Completely disassemble Polydentia Products if applicable
- We suggest to use a cleaning solution such as, but not limited to:
 - Dürr Dental ID212Or cleaning solutions based on Quaternary Ammonium compounds. Check the pH-value of the cleaning solutions; it must be as neutral as possible.
- Follow carefully the instructions for use given by the manufacturer of the selected cleaning solutions. Please note that too strong concentration or too long exposure time may damage Polydentia Products.
- Change the cleaning solutions regularly according to the instructions.
- Always rinse carefully Polydentia Products under running water after cleaning.

B. Ultrasonic bath

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Completely disassemble Polydentia Products if applicable.
- Follow carefully the instructions for use of the ultrasonic bath device.
- Service and clean the device frequently, according to the instructions for use.
- Ensure that maximum load is not exceeded.
- We recommended to use a cleaning solutions such as, but not limited to:
 - Dürr Dental ID212
 - Or cleaning solutions based on Quaternary Ammonium compounds. Check the pH-value of the cleaning solutions; it must be as neutral as possible.
- Follow carefully the instructions for use given by the manufacturer of the selected cleaning solutions. Please note that too strong concentration or too long exposure time may damage Polydentia Products.
- Change the cleaning solutions regularly according to the instructions.
- Remove the devices from the ultrasonic bath immediately after cleaning.
- Always rinse carefully Polydentia Products under running water after cleaning.

C. Thermal Disinfector

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Completely disassemble Polydentia Products if applicable.
- Follow carefully the instructions for use given by the manufacturer of the thermal disinfector and of the selected program; and prefer cleaning agents with corrosion protector.
- Inspect the device regularly according to the instructions for use. Service and clean the device frequently, according to the instructions for use.
- Ensure that maximum load is not exceeded.
- Ensure that the selected washing programme and cleaning agents are effective enough.
- Remove Polydentia Products from the thermal disinfector immediately after cleaning.

5. Drying

Before sterilizing in the autoclave, Polydentia Products must be thoroughly dried. The normal tap water may contain different chemical ingredients which may cause chemical reactions when autoclaving.

6. Checking

After cleaning, check that Polydentia Products are clean by rubbing the surface with an eraser or finger.

WARNING! It is extremely important to check that the products are truly clean and dry before the sterilizing procedure.

7. Sterilization

Sterilize Polydentia Products **only** in a steam autoclave using distilled water and at the temperature specified on the relevant symbol. **Do not use chemical or cold sterilization.**

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Follow carefully the instructions for use of the autoclave.
- Inspect the device regularly according to the instructions for use. Service and clean the device frequently, according to the instructions for use .
- Ensure that maximum load is not exceeded.
- Place all the medical device components in an autoclave tray, or in a sterilization disposable pouch (recommended).
- Plastic parts must not touch the autoclave walls, since the temperature there might be higher.
- We recommend to use always the following cycle (Prion Cycle):
Temperature of sterilisation: 134°C
Pressure of sterilisation: 2 bar
Time of sterilisation: 20.5 min
Drying time: 20 min
- Remove Polydentia Products from the autoclave immediately after sterilisation.

8. Storage

- Store sterilized product in a dry and dust free location.
- Keep product sealed in the sterilization pouch until ready to use.

DISCLAIMER:

This high quality product is manufactured by Polydentia SA, in Mezzovico, Switzerland. The content of this packaging has been inspected, tested and approved by our quality officers.

Materials used during our manufacturing process are designed to provide durability and performance. Polydentia constantly looks for the highest quality plastics for its product; nevertheless please take into account that plastics may alter over time. Sterilization or exposure to chemicals can accelerate this deterioration.

For the longest lifetime of your premium quality product always clean and sterilize it following these instructions.

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