Instruction for use

myQuickmat Forte kit

sectional
matrix
systems

Medical Devices Kit for professional use only - REV 01 06.2021

instruction for use

A - Description

This kit features myRing Forte, a sectional matrix ring that adapts to both teeth of different crown heights and wide cavity preparations, thanks to its extremities (myTines Small, myTines Medium and myTines Large). In case of dissimilar tooth morphologies, the autoclavable and replaceable myTines can be used interchangeably to enhance matrix adaptation. The kit comprises a selection of myWedge, high-quality, anatomical interdental plastic wedges. Their hollow V-shape allows them to compress and adapt to the interproximal anatomy of the teeth, preventing damage to the papilla. LumiContrast sectional matrices complete the kit's matrix assortment. These dark-blue matrices significantly increase contrast and reduce glare effect, especially while using loupes or microscopes.

B - Composition

myRing Forte: Stainless steel
myTines: thermoplastic
LumiContrast matrices: Stainless steel
myWedge: thermoplastic
myForceps: stainless steel

C - Indications

myQuickmat Forte kit is indicated for strong tooth separation in order to create contact points in the restoration of Class II cavities. myRing Forte's anatomically shaped, reusable tines ensure perfect matrix adaptation and myWedge forms the matrix band to the cervical area of the tooth

D - Contraindications

None Known.

E - Side effects

There are no known side effects. Promptly report any serious incident by contacting your distributor or Polydentia.

F - Interactions with other medical devices or products

There are no known interactions.

G - Information for the patient

No information for the patient related to myQuickmat Forte kit is necessary.

H - Preservation and storage

Store in a clean and dry environment.

L - Warnings

Due to the presence of small and sharp components, the use of a rubber dam is mandatory to avoid accidental ingestion or aspiration. Always ensure myRing Forte using dental floss. In case of accidental inspiration of a myQuickmat Forte kit component (matrix, ring and its extremities, wedge) seek immediate medical attention. In case of ingestion of a myQuickmat Forte kit component, it is strongly recommended seeking for medical assistance and locate the ingested components. If illness develops after ingestion, consult a doctor immediately. All Polydentia reusable Products are supplied non-sterile and must be appropriately disinfected, cleaned and sterilized before the first use. Please refer to the Cleaning, disinfection and Sterilization Instructions at www.polydentia.com.

M - Application

See the following pictograms.

N - Pack sizes available

See "Order now" section.

O - Disposal

Because of their composition, Matrices and wedges can not be sterilized and are therefore single use devices. Do not reuse them. Properly dispose devices contaminated with potentially infectious substances of human origin, in accordance with local regulations.

P - Ring reconditioning

For a proper, effective and safe reconditioning of the myRing Forte and their extremities read the Cleaning, Disinfection and Sterilization Instructions at www.polydentia.com.

QR instructions
Cleaning and
Sterilization
Guide or go to
www.polydentia.ch

application

1. Prepare the cavity

Remove the carious lesion and prepare the tooth with a conservative Class II preparation. Even if the contact was not totally broken the matrix can be placed.

2. Place the matrix

Select an appropriate matrix band according to the tooth morphology and restoration type. Place the contoured sectional matrix with the concave edge towards the occlusal margin using your fingers or with tweezers.


3. Wedging

Select the appropriate sized mywedge to close the gingival margin, and insert it firmly. When needed maintain a finger on the matrix band to avoid dislodging during wedging.


4. Place myRing Forte

Expand myRing Forte using Polydentia's myForceps; enlarge it just enough to position it on the tooth. Carefully place myRing Forte over the wedge as far gingivally as possible.


5. Burnish the matrix

For an excellent contact, burnish the band against the adjacent tooth using a ball burnisher or tweezers.


6. Restore the cavity

Restore the cavity according to the filling material manufacturer's instructions.

7. Remove the sectional matrix system

Remove myRing Forte, the wedge and the matrix. Our suggestion for an easy matrix removal is to use the special Polydentia myQuickmatrix Forceps.

8. Finishing and polishing



With thin burs remove any small excess of material then polish the restoration according to the finishing and polishing instruments manufacturer's recommendations.

tips & tricks using myQuickmat Forte kit



Tines removal

Tines must be disassembled before cleaning and sterilizing. To remove the tines correctly, it is important to preserve the internal anatomical wings. Make sure to pull the tines by grasping the external part. For more information about cleaning and sterilizing refer to the document "cleaning and sterilization procedure"

Tines placement on myRing Forte


Tines are replaceable and supplied in pairs. You can easily detach the tines, avoid finger pressure on the internal anatomical wings to preserve them. Make sure tines are placed on myRingForte as shown in the picture below.

MOD and multiple tooth restorations

MyQuickmat Forte rings are designed to be stacked over one another to perform a MOD and multiple tooth restorations.


1. Prepare the cavities as described in the instructions for use, step 1.
2. Place matrices and wedges on both mesial and distal surfaces of the tooth being restored
3. Position myRing Forte
4. Complete the restoration as described in the instructions for use, steps 5 to 8.

Restoration of adjacent cavities

MyRing Forte tines are shaped to perfectly adapt the mesial and distal surfaces of two adjacent teeth. This feature allows the restoration of two adjacent cavities at the same time.


1. Prepare the cavities as described in the instructions for use, step 1.
2. Place matrices and wedges on both teeth being restored
3. Position myRing Forte
4. Complete your restoration as described in the instructions for use, steps 5 to 8.

myRing Forte Reset

If myRing Forte is over-expanded it can be retensioned by using Polydentia's myForceps as follows:

- remove the tines first
- place myRing Forte within the notches on the forceps
- squeeze forceps closed to return myRing Forte in its perfect shape.

symbols

- REF** Item number
- LOT** Batch code
-  Use by
-  Consult instructions for use
-  Do not reuse
- R_x Only** Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner
-  Manufacturer
- CE** CE mark
-  Ultrasonic bath
-  Washer-disinfector for thermal disinfection
-  Sterilizable in a steam sterilizer (autoclave) at the temperature specified
- MD** Medical Device
- EC REP** European Authorised Representative

reprocessing of Polydentia reusable products

Disinfection,
Cleaning &
Sterilisation
instructions

polydentia
swiss manufacture 

Index

EN - English	2
DE - Deutsch	6
FR - Français	11
IT - Italiano	16
ES - Español	21

Reprocessing of Polydentia reusable Products

Disinfection, Cleaning & Sterilisation instructions

Basic Principles:

- Polydentia Products are supplied non-sterile and must be appropriately disinfected, cleaned and sterilized before the first use.
- Polydentia Products must to be disinfected, cleaned and sterilized prior to each use.
- Effective cleaning is an indispensable requirement for proper Polydentia Products sterilization.
- Consider the legal provisions valid for your country as well as to the hygienic instructions of the Dental Practice.
- The user is responsible for the sterility of Polydentia Products. Therefore, please ensure that only validated procedures are used for disinfection, cleaning and sterilization. The sterilization equipment must also be maintained and checked regularly, as well as the validated parameters applied to each disinfection, cleaning and sterilization cycle.
- Materials may alter over time. Sterilization or exposure to chemicals can accelerate this deterioration. Always check your Polydentia Product for wear and damage before use.
- Avoid contact of dissimilar metals at any point during the disinfection, cleaning and sterilization process.

Protection of Staff Members:

Used and contaminated Polydentia Products must be handled with appropriate Personal Protective Equipment.

1. Pre-treatment

Before processing Polydentia Products remove coarse impurities on the products immediately after application.

Do not use a metal brush. A metal brush will cause scratches or damages on the products surface.

- ▶ **In case of using a thermal disinfector the steps in this guidance (see below) “2. Disinfection” and “3. Rinsing” are unnecessary, provided that the dirty Polydentia Products have not been stored dry longer than 2 hours.**

2. Disinfection

- Do not leave dirty Polydentia Products dry – put them into a disinfecting liquid as soon as possible.
- We suggest to use an enzymatic pre-cleaner disinfectant such as, but not limited to:
 - **Dürr Dental ID215**
 - Or a pre-cleaning disinfectants based on Quaternary Ammonium compounds.
- Follow carefully the instructions for use given by the manufacturer of the selected disinfectant. Please note that too strong concentration of the disinfectant or too long exposure time may damage the products.
- Change the disinfectant regularly according to the instructions.

3. Rinsing

- Immediately after the disinfection, Polydentia Products must be rinsed carefully under running lukewarm water.

4. Cleaning

Always clean Polydentia Products, disinfection and rinsing are not sufficient.

Cleaning alternatives:

Manual (A), Ultrasonic bath (B) Thermal Disinfector (C)

If possible, an automatic procedure in a thermal disinfector should be used for cleaning of Polydentia Products. A manual procedure, even in case of application of an ultrasonic bath, should only be used if an automatic procedure is not available or if such a method is not compatible with specific materials; in this case, the significantly lower efficiency of a manual procedure must be considered.

A. Manual Cleaning

- **Completely disassemble Polydentia Products if applicable**
- We suggest to use a cleaning solution such as, but not limited to:
 - **Dürr Dental ID212**
 - Or cleaning solutions based on Quaternary Ammonium compounds. Check the pH-value of the cleaning solutions; it must be as neutral as possible.
- Follow carefully the instructions for use given by the manufacturer of the selected cleaning solutions. Please note that too strong concentration or too long exposure time may damage Polydentia Products.
- Change the cleaning solutions regularly according to the instructions.
- Always rinse carefully Polydentia Products under running water after cleaning.

B. Ultrasonic bath

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Completely disassemble Polydentia Products if applicable.
- Follow carefully the instructions for use of the ultrasonic bath device.
- Service and clean the device frequently, according to the instructions for use.
- Ensure that maximum load is not exceeded.
- We recommended to use a cleaning solutions such as, but not limited to:
 - Dürr Dental ID212
 - Or cleaning solutions based on Quaternary Ammonium compounds. Check the pH-value of the cleaning solutions; it must be as neutral as possible.
- Follow carefully the instructions for use given by the manufacturer of the selected cleaning solutions. Please note that too strong concentration or too long exposure time may damage Polydentia Products.
- Change the cleaning solutions regularly according to the instructions.
- Remove the devices from the ultrasonic bath immediately after cleaning.
- Always rinse carefully Polydentia Products under running water after cleaning.

C. Thermal Disinfector

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Completely disassemble Polydentia Products if applicable.
- Follow carefully the instructions for use given by the manufacturer of the thermal disinfector and of the selected program; and prefer cleaning agents with corrosion protector.
- Inspect the device regularly according to the instructions for use. Service and clean the device frequently, according to the instructions for use.
- Ensure that maximum load is not exceeded.
- Ensure that the selected washing programme and cleaning agents are effective enough.
- Remove Polydentia Products from the thermal disinfector immediately after cleaning.

5. Drying

Before sterilizing in the autoclave, Polydentia Products must be thoroughly dried. The normal tap water may contain different chemical ingredients which may cause chemical reactions when autoclaving.

6. Checking

After cleaning, check that Polydentia Products are clean by rubbing the surface with an eraser or finger.

WARNING! It is extremely important to check that the products are truly clean and dry before the sterilizing procedure.

7. Sterilization

Sterilize Polydentia Products **only** in a steam autoclave using distilled water and at the temperature specified on the relevant symbol. **Do not use chemical or cold sterilization.**

- **Please assure that Polydentia Products have the relevant symbol on the packaging.**
- Follow carefully the instructions for use of the autoclave.
- Inspect the device regularly according to the instructions for use. Service and clean the device frequently, according to the instructions for use .
- Ensure that maximum load is not exceeded.

- Place all the medical device components in an autoclave tray, or in a sterilization disposable pouch (recommendet).
- Plastic parts must not touch the autoclave walls, since the temperature there might be higher.
- We recommended to use always the following cycle (Prion Cycle):
Temperature of sterilisation: 134°C
Pressure of sterilisation: 2 bar
Time of sterilisation: 20.5 min
Drying time: 20 min

- Remove Polydentia Products from the autoclave immediately after sterilisation.

8. Storage

- Store sterilized product in a dry and dust free location.
- Keep product sealed in the sterilization pouch until ready to use.

DISCLAIMER:

This high quality product is manufactured by Polydentia SA, in Mezzovico, Switzerland. The content of this packaging has been inspected, tested and approved by our quality officers. Materials used during our manufacturing process are designed to provide durability and performance. Polydentia constantly looks for the highest quality plastics for its product; nevertheless please take into account that plastics may alter over time. Sterilization or exposure to chemicals can accelerate this deterioration. For the longest lifetime of your premium quality product always clean and sterilize it following these instructions.

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